


# FREEPRINT®

## »» SPLINTMASTER «« TAFF & FLEX

Ultimate solution for splints, repositioners, mouth- & nightguards

Now  
available!



- ✓ Extreme impact strength
- ✓ Flexible and robust
- ✓ Oral resistant
- ✓ Long-term colour stable
- ✓ Odour- & tasteless



# FREEPRINT® SPLINTMASTER

Light-curing formulation for 3D printing of flexible high-end splints. Available in two levels of flexibility: Taff & Flex.

The biocompatible material is characterized by extreme impact strength. Splints made of FREEPRINT® SPLINTMASTER are long-term stable and robust, without becoming brittle. Thanks to its high elastic recovery, the splints, always return to their original shape. The perfectly adjusted properties guarantee highest, tension free wearing comfort and easy insertion and removal. The clear-transparent flexible splints are odour- and tasteless and free from MMA, THF-MA and TPO.

FREEPRINT® SPLINTMASTER meet all requirements for a **Class IIa Medical Device according to the new MDR** and is currently undergoing **FDA 510(k) clearance**.

- ✓ Tension free wearing comfort
- ✓ Effortless insertion & removal
- ✓ Easy to clean
- ✓ Clear-transparent
- ✓ MMA-, THF-MA- & TPO-free



REF 04433  
REF 04432

FREEPRINT® SPLINTMASTER TAFF, 1.000g  
FREEPRINT® SPLINTMASTER FLEX, 1.000g

