



PREMIUM QUALITY



MADE IN GERMANY

5TH GENERATION
ZIRCONIA



NEW!

CopraHyperion MD PRO

Discover yourself the

Smart Gradient Technology



without
loss of
translucency



Super-Speedsintering**

60 min* (up to 3 units)

** Super-Speedsintering under vacuum
in the Ivoclar© Programat CS6

Speedsintering*

2 hrs. (up to 3 units)

Our masterpiece - innovation meets aesthetics

The *Smart Gradient Technology* enables a stepless color and translucent gradient. The new CopraHyperion MD PRO is a further development of our successful CopraSupreme Hyperion product range.

Isostatically pressed with a lifetime warranty.

* can lead to reduced bending strength

	CopraHyperion MD PRO
Zirconia	5th generation, 100% Tosoh
Isostatically pressed	yes
Bending strength	~ 900-1200 MPa
Class	IIa MDR
Aesthetics	layerless transition multilayer
Translucency	~ 45% - 49%
Veneering	all porcelains suitable for zirconia
Colours	A1-D4 and 0M1- 0M3
Heights	12mm, 14mm, 16mm, 18mm, 20mm, 25mm, 30mm in Ø 98 and Ø 95

Areas of
application



veneer



inlay



anterior crown



posterior crown



full arch
crown bridge



full crown
bridge



screw-retained
bridge