

GumQuick Pro



Technical Data Sheet

Product description: Universal gingival mask material for C&B technique and implantology

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Vinyl-containing polydimethylsiloxane
2	Polymethylhydrogensiloxane
3	Silicic acid
4	Filling materials
5	Pigments
6	Platinum catalyst

Properties:

Color: pink opaque

Processing time: $\geq 1 \text{ min } 30 \text{ s} \pm 10 \text{ s}$ at 23 °C

Setting time: 5 min $\pm$ 30 s at 23 °C

Final hardness: 72 $\pm$ 3 Shore A (DIN ISO 48-4)

Consistency: light-bodied - Type 3 (DIN EN ISO 4823)

Elastic recovery: >99.5 % (DIN EN ISO 4823)

Linear dimensional change: <0.5 % (DIN EN ISO 4823)

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@dreve.de

dentamid.dreve.de